



# Devnote-Cx43 liposomes

⚙ Status	Shipped
≡ Tags	
👤 Researchers	Ⓜ Yen-Yu Hsu
≡ Keywords	
✚ Summary	The research investigates the localization of Cx43-EGFP in liposome membranes and its functionality as a hemichannel using dye leakage assays. Initial results show successful expression and localization of Cx43-EGFP, although expression levels are low. Observations indicate that liposomes with Cx43 exhibit reduced fluorescence and yield compared to controls, suggesting potential dye leakage through Cx43 channels.



## Overview

*Summary of the research or experiment, particularly the **why** of the work.*

We're trying to see

(1) If PURE expressed Cx43 can localize to the liposome membranes

(2) If the Cx43-EGFP can function as a hemichannel using dye leakage assay



## Contents

@August 27, 2025

@August 29, 2026



## Notebook

### MATERIALS

#### Lipids

POPC: 16:0-18:1 PC (POPC)

Chol: cholesterol (plant)

Deuterated lipid: 16:0-d31-18:1 PC

Rhod-PE

Glucose

**@August 27, 2025**

We want to see if Cx43-EGFP can be expressed using NEB PURE and localize to liposome membranes.

### REACTION SETUP

**Lipid-in-oil solution (~5mg/ ml):** 74.5% POPC + 29.95% Chol + 0.05% Rhod



1. Add 1 ml mineral oil in the 2mL small glass jar
2. Add lipids shown in the above table into the glass jar on top of the mineral oil

3. Vortex the lipid-in-oil mixture for 10 secs
4. Put the glass jar in the bead-loaded hot bath at ~55c for 3 hrs (keep the jar uncovered without lids)
5. Place the jar (with lid) containing lipid-in-oil solution at RT for 10 mins before using


**Outer solution:** 10 mM HEPES-KOH and 10 mM NaCl solution in 700 mM Glucose

**Inner solution:**


Sample	+ Cx43	- Cx43
Sol A	8	8
Sol B	6	6
RNase Inhibitor	1	1
DNA (Cx43-deGFP)	2	0
Optiprep	0.67	0.67
H2O	2.33	4.33
<b>Total</b>	20	20

## RESULTS

### ▼ Control




561: 80 ms




488: 200 ms



561+488



561



488



561+488



561



488



561+488



561




488




561+488


▼ Sample




561




488




561+488



561




488



561+488




3 hrs of expression



3 hrs of expression


## ▼ Cells - Samples





561 channel: exp 80 ms; 488 channel: exp 200 ms.

## ▼ Data



## CONCLUSION

Cx43-EGFP was successfully expressed using the NEB PURE system and localized to liposome membranes. However, the expression level appears to be low, as indicated by dim green fluorescence. In addition, some liposomes didn't show green fluorescent ring.

## OBSERVATIONS

Although the overall expression level was low, we observed clear localization of Cx43-EGFP on liposome membranes. In contrast, control liposomes lacking Cx43

showed no green fluorescent rings under identical imaging conditions using the 488 nm channel (200 ms exposure).

**@August 29, 2026**

Repeat the experiment from @August 28, 2025

#### REACTION SETUP

**Lipid-in-oil solution (~5mg/ ml):** 74.5% POPC + 29.95% Chol + 0.05% Rhod



1. Add 1 ml mineral oil in the 2mL small glass jar

2. Add lipids shown in the above table into the glass jar on top of the mineral oil
3. Vortex the lipid-in-oil mixture for 10 secs
4. Put the glass jar in the bead-loaded hot bath at ~55c for 3 hrs (keep the jar uncovered without lids)
5. Place the jar (with lid) containing lipid-in-oil solution at RT for 10 mins before using

**Outer solution:** 10 mM HEPES-KOH and 10 mM NaCl solution in 700 mM Glucose

**Inner solution:**

Sample	+ Cx43	- Cx43
Sol A	8	8
Sol B	6	6
RNase Inhibitor	1	1
DNA (Cx43)	2	0
Optiprep	0.67	0.67
Alexa Fluor 647 (1mM)	0.1	0.1
H2O	2.23	4.23
<b>Total</b>	<b>20</b>	<b>20</b>

## RESULTS

### ▼ Figures



Liposome fluorescence intensity



Top 1% liposomes

▼ Sample



Startpoint

Endpoint





Timeseries: Every 10 mins

▼ Control



Startpoint

Endpoint



Timeseries: Every 10 mins



Timeseries: Every 10 mins

## OBSERVATIONS

- Liposomes containing Cx43 show lower yield and smaller average size compared to liposomes without Cx43.
- The average dye fluorescence of liposomes with Cx43 is markedly reduced relative to controls, suggesting dye leakage through Cx43 channels.
- A higher proportion of non-fluorescent liposomes is observed in the Cx43-containing samples compared to those without Cx43.